pempa®

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PL EN

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
USER MANUAL


COMPONENTS

1. Compressor housing
2. ON/OFF switch
3. Power cord
4. Air outlet
5. Nebulizer holder
6. Air-vent openings


SET CONTAINS:

9. Air tube (PVC)
10. Nebulizer (PP)
11. Mouthpiece (PP)
12. Nosepiece (PP)
13. Adult mask (PVC)
14. Pediatric mask (PVC)
15. Spare air filters
16. Accessories case

INSTRUCTIONS FOR USE

- The device must be checked before each use, in order to detect possible functioning anomalies and/or damages due to transport and/or storage. During inhalation, sit upright and relaxed at a table and not in an armchair, in order to avoid compressing your respiratory airways and impairing the treatment effectiveness.
- The accessories must be used only with a single patient, it is not recommended to use them with several patients.
1. After unpacking the device, check it for visible damages or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check for accessories integrity.
 2. Before using the device, proceed with the cleaning operations as described in the "CLEANING AND MAINTENANCE" section.
 3. Open the nebulizer by turning counterclockwise the top (**Pic. A**).
 4. Make sure that the medicine conduction cone is properly fitted on the air conduction cone inside the nebulizer (**Pic. B**).
 5. Put the prescribed quantity of medicine into the nebulizer (**Pic. C**).
 6. Close the nebulizer by turning clockwise the two parts, paying attention that they are thoroughly sealed (**Pic. D**).
 7. Connect one end of the air tube to the nebulizer (**Pic. E**) and the other end to the air outlet (**Pic. F**) on the device.
 8. If you need to use the mask, insert it directly onto the nebulizer (**Pic. G**).
 9. Plug the device into the wall socket, making sure that the main supply corresponds to the electrical rating of the device.
 10. To start the treatment, set the I/O switch into the "I" position.
 11. **ATTENTION: this device is for intermittent use: 30 min ON / 30 min OFF.**
 11. Inhale the aerosol solution using the prescribed accessory.
 12. When the treatment has been completed, switch the unit off y setting the I/O into the "O" position and disconnect the plug from the wall socket.
 13. Wash the nebulizer and its accessories as described in the "CLEANING AND DISINFECTION OF THE ACCESSORIES" section.
 14. The device requires no calibration.
 15. No modification to the device is permitted.

IMPORTANT: READ CAREFULLY THE INSTRUCTIONS FOR USE BEFORE OPERATION.

- INTENDED USE**
- PEMPA NEB200 is an aerosoltherapy system suitable for domestic use.
- PEMPA NEB200 is designed for the production of compressed air to operate a nebulizer kit for the production of medical aerosol for respiratory disorders. Patient population: The device is intended for use with children from 2 years old, adolescent and adult patients. Intended users: the use of the device does not require a specific knowledge or professional ability. The patient is the intended operator except in case of child and patient that required special assistance
- WARNINGS:**
1. Use the device only as described in this manual and therefore as an aerosoltherapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
 2. Keep this manual for future reference.
 3. Do not operate the unit in presence of any anaesthetic mixture inflammable with oxygen or nitrogen protoxide.
 4. The correct functioning of the equipment can be affected by electromagnetic interferences which exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it to a different power socket.
 5. In case of failure and/or malfunction, read the "POSSIBLE PROBLEMS AND HOW TO SOLVE THEM" section. Do not handle or open the compressor housing.
 6. For repair operations address only to a technical service centre authorized by the manufacturer and require the use of original spare parts. The non-observation of the above mentioned indications can compromise the device safety.
 7. Comply with the safety regulations concerning the electrical devices and in particular:- use only original accessories and components;
 - never submerge the unit in water;
 - never wet the device, it is not protected against water penetration;
 - never touch the unit with wet or moist hands;
 - do not leave the unit exposed to the weather elements;
 - place the unit on a stable and horizontal surface during its operation;
 - the use of this device by children and disabled requires always the close supervision by an adult with full mental faculties;
 - do not pull the power cord or the device itself to unplug it from the power socket;
 - the power plug is the separation element from the grid power; keep the plug accessible when the device is in use.
 8. Before plugging in the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
 9. In case the power plug provided with the device does not fit your wall socket, address to qualified personnel for the replacement of the plug with a suitable one. In general, the use of adapters, simple or multiple, and/or extension cables is not recommended. If their use is indispensable, it is necessary to use types complying with safety regulations, paying however attention that they do not exceed the maximum power limits, indicated on adapters and extension cables.
 10. Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not operated.
 11. The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to persons, animals or things, for which the manufacturer cannot be held responsible.
 12. The power cord of this device cannot be replaced by its use. In case of a power cord damage, address to a technical service centre authorized by the manufacturer for its replacement.
 13. The power supply cord should always be fully unwound in order to prevent dangerous overheating.
 14. Before performing any maintenance or cleaning operation, turn off he device and disconnect the plug from the main supply.
 15. Some parts of the unit are so small that they may be swallowed by children; keep the equipment out from children's reach.

16. If you decide to use the device any longer, it is recommended to dispose of it according to the current regulations.
17. Make sure to:
 - use this device only with medicines prescribed by your doctor;
 - make the treatment using only the accessory recommended by your doctor depending on the pathology ;
 - use the nosepiece accessory only if expressly indicated by your doctor and paying attention to NEVER introduce the bifurcations in the nose, but only bringing them as close as possible.
18. Check in the medicine package leaflet for possible contraindications for use with common aerosolotherapy systems.
19. To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.
20. Do not position the equipment so that it is difficult to operate the disconnection device.
- 21.Nebulizer and accessories are single patient reuse. Device is multi-patient reuse
22. Nebulizing system is not suitable for use in anaesthetic breathing system or a ventilator breathing system
23. Device is designed to nebulize solution and suspension liquids

CLEANING AND MAINTENANCE

- The cleaning of the device must be carried out by using a soft and dry cloth and non-abrasive cleansers.
- ATTENTION:** during cleaning operations, make sure the internal parts of the unit are not in contact with liquids and that the power plug is disconnected.
- CLEANING AND DISINFECTION OF THE ACCESSORIES**
- Follow carefully the cleaning and disinfection instructions of the accessories as they are very important for the device performances and the therapy success.
- 1. First use and after each treatment:**
- Disassemble the nebulizer by turning counterclockwise the top and remove the medicine conduction cone. Wash the components of the disassembled nebulizer, the mouthpiece, the masks and the nosepiece by using tap water; dip in boiling water for 5 minutes. Reassemble the nebulizer components and connect it to the air-outlet (Pic.G), switch the device on and let it work for 10-15 minutes.
- ATTENTION:** The masks and air tube must be wash ed with warm water.
- ATTENTION:** Do not boil nor autoclave the air tube and masks.
- It's possible to use cold disinfecting liquids following the manufacturer's instructions.

REPLACEMENT OF THE NEBULIZER

- The nebulizer must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebulizer nozzle is obstructed by dry medicine, dust, etc. it is recommended to replace the nebulizer after a period between 6 months and 1 year depending on the usage.
- Use original nebulizers only.**
- REPLACEMENT OF THE AIR FILTER**
- In normal usage conditions, the air filter must be replaced approximately after 200 working hours or after each year.
- It is recommended to periodically check the air filter (10 – 12 treatments) and, if the filter shows a grey or brown colour or is wet, replace it. Extract the filter (Pic. A-5) and replace it with a new one. Do not try to clean the filter for reusing it. The air filter shall not be serviced or maintained while in use with a patient.
- ATTENTION:** Use original filters only. Do not use the device without filter.
- POSSIBLE PROBLEMS AND HOW TO SOLVE THEM**
- Contact an authorized Customer Service Centre in case you need help about the use and the maintenance of the equipment.
- PROBLEMS:**
- 1.The device does not switch on2.The device does not nebulize or nebulizes weakly
- SOLUTIONS:**
- 1.Make sure the power plug is firmly fitted to the wall socket
 2. Make sure that the device has been operating within operating limits indicated in this manual (30 min ON / 30 min OFF).
 - 3.Make sure that the ends of the air tube are tightly fitted to the main unit and nebulizer;
 - 4.Verify that the nebulizer is not empty or has been filled with the poper quantity of medicine (max 8 ml).
- Verify that the nebulizer nozzle is not obstructed.
- Should the device not start working properly again, address to customer service of authorized retailer
- MAINTENANCE AND REPAIRS** In case of failure, address to qualified personnel authorized by manufacturer. Do not open the device in any case. The unit has no user-serviceable parts within and does not need internal maintenance or lubrication.

Please report any serious incident that has occurred in relation to the device injury or adverse event to the local competent authority and to the Manufacturer or to the European Authorised Representative (EC REP)

Vigilance contact point: <https://ec.europa.eu/growth/sectors/medical-devices/contacts/>

TECHNICAL SPECIFICATIONS

- Model/Type: PEMPA NEB200/GCEB30
- Power supply: 230V 50 Hz AC
- Current: ≤ 1000mA
- Minimum filling volume: 2 ml
- Maximum filling volume: 8 ml
- Nebulization rate: average 0,40 ml/min (range 0.23-0.60ml/min)
- Maximum air pressure: average 2,2 bar (range 2.0 - 3.2 bar)
- Weight: 1.5 kg
- Size: 160x92x140mm
- Noise level: average 52 dBA (max level 58dBA)
- Expected service life of device: 1000 hours
- IP protection: IP21
- Class II** device as regards protection against electric shocks.
- Nebulizer, mouthpiece and masks are type BF applied parts.
- Device not suitable for use in presence of anaesthetic mixtures inflammable with air, oxygen or nitrogen protoxide.
- Device for intermittent use (30 min ON / 30 min OFF)
- Device not suitable for use in anaesthesia or lung ventilation systems.
- NOTE: The technical specifications may change without prior notice.
- Applied standards**
- Electric Safety Standards EN 60601-1. Electromagnetic Compatibility according to EN 60601-1-2.
- Device complies with the requirements of the Medical Device Directive 93/42/EEC.

AEROSOL PERFORMANCES ACCORDING TO EN13544-1:2009

- based on adult ventilatory pattern with Sodium Fluoride (NaF)
- Aerosol output: 0.174 ml
- Aerosol output rate: 0.082 ml/min
- Percentage of fill volume emitted per min: 4%
- Residual Volume: 1.576 ml

- Particle size (MMAD): 3.16 µm
- GSD (geometric standard deviation): 2.27
- RF (respirable fraction <5µm): 71.2%
- Large particle range (>5µm): 28,8%
- Middle particle range (2 to 5µm): 40,9%
- Small particle range (<2µm): 30,3%

ENVIRONMENT CONDITIONS

- Storage / Transport
- Temperature: MIN -25 °C – MAX +70 °C
- Humidity: MIN 10% RH – MAX 95% RH
- Atmospheric pressure: 700 hPa – 1060 hPa

OPERATION

- Temperature: MIN +10 °C - MAX +40 °C
- Humidity: MIN 10% RH - MAX 95% RH
- Atmospheric pressure: 700 hPa – 1060 hPa

GCEB30-IM-MDD-2020-12-14-Rev00

Tab. 1) Guidance and manufacturer's declaration - Electromagnetic emissions			
The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment			
Emission test	Compliance	Electromagnetic environment - guidance	
Emissions RF CISPR 11	Group 1	The Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class (B)	The Device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Class A		
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies		

Tab. 2) Guidance and manufacturer's declaration - Electromagnetic immunity			
The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±18 kV contact ±15 kV air	±18 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	Power supply line: ±2 kV input/output lines: ±1 kV 100 kHz repetition frequency	Power supply line: ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV line(s) to earth: ±2 kV	line(s) to line(s): ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5% UT (≥95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle ≤5% UT (≥95% dip in UT) for 5s	5% UT (≥95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle ≤5% UT (≥95% dip in UT) for 5s	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	20 A/m 50Hz/60Hz	20 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE 1 In the a.c. mains voltage prior to application of the test level.

Tab. 3) Guidance and manufacturer's declaration - electromagnetic immunity			
The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	150kHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Device including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances:
Radiated RF IEC 61000-4-3	10V/m	10V/m	80 MHz to 800MHz d=1.2√P 800 MHz to 2.7 GHz d=2.3√P

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. In cases the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Tab. 4) Recommended separation distances between portable and mobile RF communications equipment and the Model Smart			
The Model Smart is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz (out GSM and amateur radio bands) d=1.2√p	150 kHz to 80 MHz (in ISM and amateur radio bands) d=2√p	800MHz to 2.7GHz d=2.3√p
0,01	0,12	0,2	0,23
0,1	0,38	0,632	0,73
1	1,2	2	2,3
10	3,8	6,32	7,3
100	12	20	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Electromagnetic immunity

The Device is intended for use in the electromagnetic environment specified below.

The customer or the user of the device should assure that it is used in such an environment.

	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	Immunity Test Level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 - 390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380 - 390	GMR5 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 - 787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	780						
	810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA IS5, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
	930						
	1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5240						
	5785						

Notes: If necessary to achieve the **Immunity Test Level**, the distance between the transmitting antenna and the **ME Equipment** or **ME System** may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50% duty cycle square wave signal.
- c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The **Manufacturer** should consider reducing the minimum separation distance, based on **Risk Management**, and using higher **Immunity Test Levels** that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher **Immunity Test Levels** shall be calculated using the following equation:

$$E = \sqrt{P} \cdot d$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the **Immunity Test Level** in V/m.

SYMBOLS OF THE DEVICE / DEVICE SYMBOLS

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WEEE - waste electrical and electronic equipment

This product is subject to the European Directive 2012/19 / EU on Waste Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic devices together with household waste. Please find out about local regulations for the proper disposal of electrical equipment and electronic devices. Correct disposal helps to prevent negative consequences for the environment and human health.

GLOBALCARE MEDICAL TECHNOLOGY CO., LTD.
7th Building, 39 Middle Industrial Main Road, European Industrial Zone, Xiaolan Town, 528415 Zhongshan City, Guangdong Province, PEOPLE'S REPUBLIC OF CHINA.
Phone : +86 760 22589901
<http://www.globalcare.com.hk/contact/>

EC REP

Donawa Lifescience Consulting Srl
Piazza Albania, 10, 00153 Rome Italy
<https://www.donawa.com/wli/main/contatti.index>

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